

PRODUCT COMPARISON:

Bionix® Curette Tip Styles


Our Bionix® Currettes offer a variety of tip styles to match different types of wax and physician preference.

WHAT ARE YOUR CURETTE TIP OPTIONS?

With our varying curette tip options, you'll be able to find the right tool for every cerumen situation. See what sets each tip apart and how our varying styles can help you meet specific patient care needs.



SAFE EAR CURETTES™	FlexLoop® 4mm All-purpose curette for a variety of patients	LIGHTED EAR CURETTE™
	CeraSpoon® 4mm Designed for soft wax and patients who produce excessive amounts of cerumen	
	ControLoop® 4mm Designed with a “Shapleigh” tip for pulling wax straight out of the ear canal	
	AngleLoop® 4mm Designed with an angled tip for hard-to-get cerumen in older patients	
	VersaLoop® 3mm Smooth, flexible tip designed for increased control and improved patient experience	
	InfantScoop® 2mm Designed for smaller ears or use through an otoscope	
	MicroLoop® 3mm Designed with a rigid construction to remove tough, impacted cerumen	



CORRESPONDING LIGHTED EAR CURETTE™

With our Lighted Ear Curette™, you can combine illumination and magnification in one handheld tool. These features make the cerumen easy to see—giving your patients a quick procedure they can be confident in.

