# **OtoClear® Ear Irrigation Tip**



info@sanibelsupply.com sanibelsupply.com

Bionix OtoClear® Ear Irrigation allows healthcare professionals to utilize a soft, flared, 3 streamed tip that eliminates the mess, improves patient safety, and is compatible with a variety of delivery devices to meet the personal preference of every facility.





#### **PROBLEM**

Traditional ear irrigation shoots water directly into the ear canal, which can be dangerous and run the risk of causing perforation to the eardrum.

### **SOLUTION**

OtoClear® Ear Irrigation Tips direct water toward the ear canal wall, with three divergent streams, instead of directly at the Tympanic Membrane. The tip is also designed with a soft, flexible material for patient comfort and flared design to prevent backsplash and over insertion.

### **RESULTS**

A clean and safe procedure that improves efficiency and patient experience.

"Our office is very happy with the soft (OtoClear) tips. They are easier to maneuver in the ear canal and more comfortable for the patient." - Laura Griffith, MA





OTOCLEAR TIPS ARE











INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

# **OtoClear® Ear Irrigation Tip**

### PRODUCT OFFERING



## **OtoClear SprayWash Kit**

- Manual, intermittent lavage.
- Delivers greater volume compared to syringes or other manual devices.
- Efficient design allows for two-handed procedure.
- 500ml bottle volume.



### **OtoClear Ear Irrigation Tips**

- Directs water to the ear canal wall, with three divergent streams, instead of directly at the Tympanic Membrane.
- Flared design prevents backsplash and over insertion.

OtoClear Tip US Patent # 6,706,023.

Aquabot® is registered under Lunatech®, used under permission of Lunatech®, Inc. All rights reserved.

