



Description

Degreaser gel for sensitive skin preparation for EEG, EP, ECG, Stress Test, Holter, Telemetry recordings. Guarantees low impedance and increased electrode adhesion.

- Viscous paste.
- Water based.

UMDNS 11425

Composition

Ingredients

Water, preservative, propylene glycol, acrylic polymer, fine pumice, sodium hydroxide.

Packaging

PE

Safety

Biocompatibility (UNI EN ISO 10993)

Yes

Microbiological test

Challenge test

Chemical and physical properties

pH

$6 \pm 0,4$

Viscosity

60.000 ± 10.000 CPS

Specifications

Generic code

EVERIXXXXXXXXXX

Package

Bottle 160gr

Environment

RoHS II

✓

REACH

✓